ENGLISH

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of about 1.5 to 2.8 mm for the very first time, depending on the VISCOJECT™ model and operating technique used, thanks to the new VISCOJECT™ hydraulics system. With the VISCOJECT™ injection system the majority of all foldable, one-piece acrylic lenses with an optical diameter of 6.5 millimeters can be injected, provided that the safe injection of the lens using the particular VISCOJECT™ injection system was validated in advance on the part of the lens manufacturer according to ISO 11979. The VISCOJECT™ injection system consists of the sterile VISCOJECT™ single-use injector with silicone tip (Fig. 1) and the sterile VISCOJECT™-BIO or VISCOJECT™ single-use cartridge (Fig. 2).

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APPLICATION
Folding and injecting of one-piece foldable lenses into the capsular bag or sulcus after extra-capsular cataract extraction. Three-piece foldable intraocular lenses are not suited for implantation using the VISCOJECT™ injection system. **INSTRUCTIONS FOR USE**
1. Open the blister in the sterile area and remove the cartridge. 2. Fill the cartridge tunnel of the VISCOJECT™ cartridge and the loading chamber of the cartridge sufficiently with viscoelastic solution (Fig. 3). In addition apply a drop of the viscoelastic solution directly to the silicone tip. Exclusive use of BSS solution is not allowed. **IMPORTANT:** When using hydrophobic lenses, let the visco-

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elastic solutions work for about 30 seconds to achieve sufficient hydration of the lubricating film. After 30 seconds the lens should be injected immediately. Viscoelastic materials may lose their lubricating properties when exposed to air for a longer period of time. 3. Place the foldable lens on the central hinge of the opened cartridge symmetrically as shown in Fig. 4. Align one-piece lenses in the direction of the cartridge regardless of the one-piece design. Asymmetric insertion of the lens can lead to rotation of the lens in the injector. 4. Press the lens into both guide rails with rounded sterile tweezers and gently move the lens back and forth horizontally in order to ensure that it can move freely (Fig. 5).

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5. Ensure that the lens is correctly aligned and that the cartilage sides close. Make sure that neither the edge of the lens nor the haptic are wedged in. In the case of lenses with C-Loop haptics, align these using sterile tweezers in the direction of the optics (Fig. 6). 6. Press the sides of the cartridge together until the click-lock mechanism engages (Fig. 7). 7. Insert the cartridge into the front end of the injector. Push the cartridge into the very front position (Fig. 8). 8. Carefully push the plunger forward and ensure that the silicone tip correctly enters the loading chamber. If the silicone tip cannot be introduced into the loading chamber, retract the plunger to the starting position and align the

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silicone tip straight using sterile tweezers. Continue to push the cushion until the inner spring begins to press together. Pull the plunger back a few millimetres and then push it forward again. This step ensures that the lens is correctly grasped. 9. Push the sides of the front part of the cartridge point outside of the eye. This prevents pressure generated by the VISCOJECT™ system and an increased quantity of viscoelastic solution from being transferred to the eye. 10. Guide the point of the cartridge through the incision and push it across the iris to the near edge of the pupil. 11. Slowly press the plunger forward (the lens is pressed forward almost without contact as a result of the pressure

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created by VISCOJECT™). Depending on the lens model, support the correct rection of the lens by slightly turning the injector. Slowly inject the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only push the plunger until the lens has emerged completely, even if the plunger is not yet at the limit. 12. Carefully position the lens, rotating it if necessary with the aid of using a suitable positioning hook. 13. Remove the viscoelastic material thoroughly from the eye and the lens with the standard irrigation and aspiration techniques.

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PLEASE ENGLISH:
Do not resterilize / prepare cartridge and injector. Single use only. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. **OPERATIONAL PROCEDURE**
The appropriate surgical techniques are the responsibility of the respective surgeon. He or she must assess the appropriateness of the relevant procedure based on his or her education and experience. **GUARANTEE AND LIABILITY LIMITATION**
The manufacturer guarantees that this product has been manufactured with the appropriate care and assumes no liability

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