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steriCLIN®

Hersteller

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3FCZB470xxx

Intended Purpose

The container filters are intended for single use in combination with container systems which are used as sterile barrier systems for terminally sterilised medical devices. The container filters consist of sterile barrier material (paper or non-woven) and their purpose is to enable sterilisation and to maintain sterility of the content of such containers.

Medical Device

According to the intended purpose, the product is a medical device.

Conformity for medical devices of class I is declared according to Medical Device Regulation (EU) 2017/745.

The packaging material is CE-marked.

Intended Users

Intended users is educated and trained staff working in the Central Sterile Supply Department of health care facilities.

Clinical Benefit

These sterile barrier systems are not directly applied to patients and do not have a medical indication. Prevention of infections is the clinical benefit.

Performance Characteristics

The performance characteristics do result from compliance with the product related standards for sterile barrier systems. The products comply with applicable standard ISO 11607. Imprinted chemical indicators comply with standard ISO 11140-1 (Type 1).

Intended Method of Sterilisation

The container filters are supplied in non-sterile condition and are suitable for the methods of sterilisation as follows:

STEAM

Important Instruction

Serious incidents with this container filters shall immediately be reported to the manufacturer and to the competent national authority relevant for medical devices.

Warnings

- **Do not re-use!** These container filters are intended for single use. Sterile barrier functionality may be compromised in case of multiple use.
- Sterile barrier functionality may be compromised by influence of damaging substances and chemicals.
- For sterilisation method other than intended, the proper functionality of sterile barrier cannot be guaranteed.
- Sterile barrier functionality may be compromised if a container filter with wrong dimensions is used conflicting the correct dimensions for a specific container.

Precautions

Do not use when the container filters or their packaging is damaged! In general, contact to chemicals, detergents or disinfectants shall be avoided.


Safe Disposal

Non-contaminated container filters can be disposed of together with household waste.

Storage Conditions and Stability

Recommended storing condition is at temperatures between 15°C and 25°C and at relative humidity level between 35% and 50%. **Keep dry! Keep away from sunlight and heat!**

The shelf life of container filters until the use together with a container and sterilisation does cover 5 years. Subsequently, the stability of the sterile barrier system is guaranteed for further 5 years of storing. Other limitations of storage period (e.g. recommendations of DIN 58953-8) have to be noted with respect to the storage conditions.

Preparation for Use

- The container type which is specified for a specific container filter must be identified. This shall include the verification of correct article number of the container filter.
- The correct fit of the container filter to a specific container needs to be checked: all openings of the container must be fully covered when the filter is placed into its correct position and fixed.

Use

- When inserting the container filter into the container, the instructions of use of the specific container, issued by the manufacturer of this container, must be followed strictly. Any handling instructions for the specific holder system of the container need to be respected to avoid damage to the filter or wrong positioning.
- Used container filters must be removed from containers prior to the re-processing (cleaning).
- As the container filters are for single use only and must not be re-used, the container filters need to be disposed after removal from the container.
- A new container filter must always be inserted into the container prior to sterilization. Damaged or used container filters must not be used.

Aseptic Presentation

For aseptic presentation, the process of opening of a container needs to be assessed. During aseptic presentation, the container filter must stay fixed in its position and must not be removed until the aseptic presentation is finished.

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