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
Rev.02-27/01/2022

Product information box for IFU Synergy, including product code (3520134), name, country (N.A.), and technical specifications for printing.

M.D. TECH S.R.L.
Via Fratelli Bandiera NCC,
80026 Casoria (NA)
www.md-tech.it

EN Instructions For Use
Synergy
Sodium hyaluronate
Product: Synergy is a viscoelastic solution containing cohesive sodium hyaluronate for injection during eye surgery.
Contents: Each pack contains 1.1 ml of viscoelastic solution supplied in a single-use sterile glass syringe, for intraocular use, sealed in a blister pack.

Table with 2 columns: Synergy concentration (1.0% to 2.0%) and Sodium hyaluronate (mg/ml) (10.0 to 20.0).


Description: The Synergy viscosurgical device is a clear, sterile, non-pyrogenic buffered saline solution containing highly purified and high-molecular-weight sodium hyaluronate obtained by bacterial biofermentation.
Characteristics: Synergy is a biocompatible solution with viscoelastic properties that lubricate, support and protect the tissues of the eye during ophthalmic surgery...

Indications: Synergy is indicated in eye surgery of the anterior segment, including cataract removal, insertion of intraocular lenses, corneal transplant, glaucoma and retinal detachment.
Dosage and administration: Cataract surgery - intraocular lens insertion. Before extracting the lens, Synergy is introduced into the anterior chamber to protect the corneal endothelium and maintain a deep anterior chamber.

Retinal detachment
Synergy should be slowly introduced into the vitreous cavity. With caution and by modulating the injection, Synergy can be used to separate the membranes from the retina for safe excision and release of traction.
Warnings and precautions: Any viscoelastic residue not removed at the end of the surgical procedure will be eliminated naturally via the trabecular system and Schlemm's canal.

Interactions: Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore, syringes pre-filled with Synergy should never come into contact with surgical instruments rinsed using these solutions...

Sterile opening technique
Peel back the Tyvek covering of the blister.
Remove the syringe from the blister and place it in the sterile field.
Assembly
Remove the plastic cap from the syringe tip by unscrewing it anticlockwise.

Attach the cannula to the syringe and ensure that it is securely inserted, screwing it on clockwise.
Verify the proper functioning of the device by gently depressing the plunger until Synergy comes out of the tip of the cannula.

Storage: Store Synergy between 2°C and 25°C. Do not freeze.
Symbols: Manufacturer, Do not reuse, Do not resterilize, Sterilized using steam, Consult instructions for use, Caution, Keep away from sunlight, Keep dry, Do not use if package is damaged.

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Table with 2 columns: Synergy concentration (1.0% to 2.0%) and Hyaluronate de sodium (mg/ml) (10,0 to 20,0).

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