ENGLISH

Medicel NAVIJECT™ Single-use cartridge set for injection of ONE-PIECE FOLDABLE ACRYLIC LENSES

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DESCRIPTION
 NAVIJECT™ injectors – made of high-grade titanium – are instruments for implanting acrylic and silicone foldable lenses. The lens can be injected into the eye gently and safely using an injector with a NAVIJECT™ cartridge. The NAVIJECT™ injection system consists of three parts:

1. The reusable NAVIJECT™ titanium injector (Fig. 1a, 1b)
2. The loading unit with silicone cushion (Fig. 2)
3. The sterile NAVIJECT™ single-use cartridge (Fig. 3)

The cartridge may only be used once and not resterilized/treated. Reuse or resterilization may reduce the product's performance, which can result in serious impairments to the patient's health and safety. Store cartridge set at room temperature.

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APPLICATION
 Folding and injecting of one-piece foldable lenses into the capsular bag or sulcus after extra-capsular cataract extraction. Three-piece foldable intraocular lenses are not suited for implantation using the NAVIJECT™ injection system.

INSTRUCTIONS FOR USE:

1. Clean, disinfect and sterilize the NAVIJECT™ injector (Fig. 1a, 1b) according to the instructions for use for the multiple-use injector "Treatment pursuant to EN ISO 17864" and place it in sterile condition in a sterile environment.
2. Open the blister package in the sterile area and remove loading unit.
3. Place the silicone tip on the injector head (Fig. 4). Make sure that the tip is upright.
4. Hold the loading unit with the NAVIJECT™ cartridge in such a manner that the rear cartridge flank can be guided with

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the index finger of the left hand. After 30 seconds the lenses are equipped with guiding rails. Push the lens onto the guide rails from the rear and position it in the middle of the loading chamber (Fig. 5).

Ensure that the lens is correctly aligned and be careful that when closing, neither the edge of the lens nor the haptics are wedged in. For reusable NAVIJECT™ titanium injectors, align these in the direction of the lens using spring tweezers.

5. Completely fill the cartridge tunnel of the NAVIJECT™ cartridge with viscoelastic solution and fill the loading chamber with sufficient viscoelastic solution (Fig. 6). In addition, apply a drop of viscoelastic solution to the silicone tip. The exclusive use of balanced salt solution [BSS] is not permitted.

IMPORTANT: When using hydrophobic lenses, let the viscoelastic solutions work for about 30 seconds to achieve suffi-

cient hydration of the lubricating film. After 30 seconds the lenses should be injected immediately. Viscoelastic materials may lose their lubricating properties when exposed to air for a longer period of time.

6. Press the cartridge sides together until the "click-lock" mechanism engages (Fig. 7).
7. Insert the cartridge into the front end of the injector. Push the cartridge into the very front position (Fig. 8).
8. Carefully push the plunger forwards and ensure that the silicone tip correctly enters the loading chamber. If the silicone tip cannot be introduced into the loading chamber, retract the plunger to the starting position and align the silicone tip straight using spring tweezers. Continue to push the cushion until the inner spring begins to press together or until the thread completely engages (screw-type). Pull the plunger back a few millimetres and then push it forward again. This

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