TEKLENS II

INTRAOCCULAR LENSES LENTILLES INTRAOCULAIRES INTRAOKULARLINSEN LENTI INTRAOCULARI LENTES INTRAOCULARES ΕΝΔΟΦΘΑΛΜΙΟΙ ΦΑΚΟΙ

CE 0459

P.N. 1040 Rev. F
Dec. 2015

TEKIA, INC.

17 Hammond, Ste. 414, Irvine, CA 92618 U.S.A
Tel: (949) 699-1300 • Fax: (949) 699-1302
www.tekia.com

European representative:

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Representante Europeo:

Rappresentante Europeo:

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CE 0459

DESCRIPTION

TEKIA, INC. intraocular lenses (IOL's) are optical lenses designed to replace the crystalline lens and to correct aphakia after cataract extraction. The Teklens II family clear IOLs are made of hydrophilic acrylic copolymer with a UV inhibitor (10% UV cut off @ 380 nm); The Teklens II family yellow IOLs are made of hydrophilic acrylic copolymer with a UV inhibitor and a blue light absorber. The 3 piece models have two polyethersulfone (PES) haptics. The available powers are between 0.0 to 35.0 D. The single piece models are available between 0.0 to 42.0 D. Consult website for minus powers.

INDICATIONS

The TEKIA, INC. Teklens II family Intraocular Lenses are designed to be implanted in the posterior chamber of the eye after the extraction of cataract in the following conditions: senile cataract, traumatic cataract, congenital cataract.

CONTRA-INDICATIONS

Implantation is not advisable when the IOL may aggravate an existing condition, interfere with the diagnosis or the treatment of pathology or present a risk to sight of the patient. Among these conditions are: children of less than one year old, uncontrolled glaucoma, rubeolic cataract, retinal detachment, atrophy of the iris, microphthalmia, developing chronic eye infections, endothelial corneal dystrophy, peroperative complications such as vitreous loss, hemorrhage or foreseeable post-operative complications.

WARNING

As with all surgical procedures, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Some of the potential complications of cataract surgery are: inflammations (iritidocyclitis, pupillary membrane, vitreous inflammation, CMO), infections

(endophthalmy), retinal detachment, pupillary block, iris hernia, flattening of the anterior chamber, hemorrhage (hyphema), corneal dystrophy, glaucoma, atrophy of the iris, IOL decentration, luxation and subluxation and precipitates on the IOL surface, calcification. TEKIA, INC. is not liable for any post-operative complications, which arise following the use of an IOL.

PRESENTATION - STERILIZATION

TEKIA, INC. Teklens II family Intraocular Lenses are individually packaged in a glass vial with an IOL holder in a single pouch and steam sterilized. Pouch and vial can only be opened in a sterile environment. A patient identification card and a set of labels are provided to insure IOL traceability and patient follow-up.

PRECAUTIONS FOR USE AND STORAGE

- Store at room temperature. Avoid high temperatures above 45°C.
- The IOL model, power and expiration date should be verified before opening.
- Do not use after expiration date indicated on package.
- Single use, do not re-sterilize. The IOLs are designed for single use only and cannot be re-sterilized, the sterility and specifications of the IOLs are only guaranteed for a single use.
- An IOL should only be implanted by an experienced surgeon trained in implantation techniques.
- Manipulating an IOL requires the use of the appropriate techniques and instruments. Any IOL damaged during handling should not be implanted.
- Before any implantation the surgeon should determine the power of the IOL required by the biometric data of the eye to be operated on and the characteristics of the IOL.
- TEKIA's A constant and ACD values (see outer box label) are recommended for the determination of IOL power using appropriate formula. Personalization of the surgeon constants is highly recommended.

- Do not use the IOL if the sterile pouch or glass vial has been compromised.
- Sterility of the IOL is guaranteed only as long as the individual sterile pouch and glass vial have not been opened or damaged.
- The hydrophilic material may absorb substance that is in contact with the lens optics (disinfectant, drug...)
- The use of hydrophilic acrylic IOLs should be avoided when procedures using intracameral air/gas or vitreous tamponade (silicone oils/gases) are anticipated.
- If a YAG laser posterior capsulotomy is performed take care that the optical break-through is produced behind the posterior capsule.

INSTRUCTIONS FOR USE

Open the protective sealed box and remove the pouch bag containing the glass vial and holder with IOL. Before opening the pouch make a final check of the IOL label for correct model and power. Begin at the chevron end of the sterile pouch and tear the pouch apart to release glass vial onto the sterile instrument tray. Open the glass vial screw cap and remove the gray stopper affixed to the holder containing the IOL. Use appropriate instrument to remove IOL from its holder. Examine lens closely. Rinse with sterile saline solution after lens is folded.

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