
PHACO-TIPS, WRENCHES, SLEEVES, TEST CHAMBERS AND I/A-INSTRUMENTS (REUSABLE)

ENGLISH

PHACO-TIPS, WRENCHES, SLEEVES, TEST CHAMBERS AND I/A-INSTRUMENTS (REUSABLE)

INSTRUCTIONS FOR USE

These instructions for use are to be read and used only in conjunction with the manufacturer's instruction manual. Failure to

consult the manufacturer's manual could result in improper use of this product.

Put pressure on the tip and screw the phaco tip onto the handle up to the stop position. Attention: Overtighten the tip can damage the multiple use / single use wrench as well as it can cause residue from the wrench on the surface of the tip.

TREATMENT ACCORDING TO EN ISO 17664

General principles: Before every use, all instruments must be cleaned, disinfected and sterilized; this applies particularly to the first use after delivery since all instruments are delivered unsterilized (cleaning and disinfection after removal of protective shipping packaging; sterilization). An effective cleaning and disinfection is an absolute necessity for effective sterilization. As part of your responsibility for the sterility of instruments during use, please ensure

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<p>be used; other sterilization procedures are not permitted.</p> <p>Steam sterilization: <ul style="list-style-type: none">• fractioned vacuum process (with sufficient product drying). The use of the less effective gravitation process must be guaranteed by an additional product-, sterilizer- and process-specific validation, and is the responsibility of the user (if necessary, longer sterilization times are required). • validated steam sterilizer pursuant to E DIN EN 13060 or DIN EN 285 • process validated according to DIN EN/ANSI AAMI ISO 17665 (valid commissioning and product-specific performance appraisal) • Sterilization temperature 134 °C +4/-13 °C (274 °F +6/-24 °F); with saturated steam pressure of 3.41 bar • Sterilization time (exposure time at sterilization temperature) of at least 20 minutes at 121 °C (250 °F) or 5 minutes at 134 °C (274 °F). <p>The flash sterilization process is fundamentally not permitted. Also, do not use any hot air sterilization, radiation sterilization, formalde-</p></p>
<p>PACKAGING</p> <p>We recommend sterilizing instruments in the sterilization tray and sterilization containers designed for this purpose; single-use packagings (single or double packaging) and/or other sterilization containers can be used, however, as long as these meet the following requirements: <ul style="list-style-type: none">• compliance with DIN EN ISO/ ANSI AAMI ISO 11607 parts 1+2 • suitable for steam sterilization (temperature resistance to at least 141 °C, sufficient vapor permeability) • sufficient protection of the instruments or sterilization packagings from mechanical damages • regular maintenance according to the manufacturer guidelines (sterilization container) <p>STERILIZATION</p> <p>For the sterilization, only the following sterilization procedures may</p></p>

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hyde or ethylene oxide sterilization, and also no plasma sterilization.

STORAGE

After sterilization, the instruments must be stored dry.

MATERIAL DURABILITY

When choosing the cleansing agents and disinfectants, make sure they do not contain the following ingredients:

- strongacids or lyes
- aluminum chloride
- halogens
- halogenated hydrocarbons
- furfural
- phenols
- methylene chloride
- nitrobenzene

Never clean any of the instruments, sterilization trays and sterilization containers with metal brushes or steel wool. The permitted cleansing agents may only be used individually. All instruments, sterilization trays and sterilization containers may not be exposed to temperatures higher than 138 °C (280 °F)!

REUSABILITY

The instruments can be reused as long as they are handled carefully and are undamaged and uncontaminated. The user takes

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responsibility of the user for the use of damaged or contaminated instruments. Any liability is excluded if the instruments are not used as directed.

OPERATIVE PROCEDURE

The appropriate surgical technique is the responsibility of the respective surgeon. Based on his education and experience, he is to evaluate the suitability of the respective procedure.

GUARANTEE AND LIMITATION OF LIABILITY

The manufacturer guarantees that this product was produced with appropriate care and shall assume no responsibility for incidental or consequential damages, losses or costs that should result directly or indirectly from the use of this product. Liability is solely limited to claim-related repairs that must be performed on the product which are clearly not attributed to incorrect handling.

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


	Wrenches and test chambers
	Phaco tips, sleeves and I/A-instruments

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	Batch code
	Keep dry
	Keep away from sunlight
	Consult instructions for use
	Manufacturer

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