EN	ENGLISH
----	---------

## Sodium Hyaluronate Ophthalmic Solution

### Information for Use

#### DESCRIPTION

S<sub>200</sub>-VISC is a clear, colorless, odorless, sterile, isotonic, viscoelastic solution for intraocular use. With its special viscoelastic properties, it maintains a deep anterior chamber and protects the corneal endothelium throughout surgery. The viscoelastic properties of the preparation allow tissue structures to be moved without causing trauma, for example in cases of prolapse of the iris. Diagnostic and therapeutic procedures may be carried out in the eye without any impairment of vision. S<sub>200</sub>-VISC protects the intraocular tissue throughout surgery and prevents adhesion and synchia formation when performing ocular surgery. Once surgery is finished, it must be removed completely by irrigation and/or aspiration. Any residue that has not been removed from the anterior chamber will be flushed out naturally via the trabecular meshwork system and Schlemm's canal. If it is not completely removed, intraocular pressure increases due to obstruction of the trabecular system.

#### BOX CONTENTS

- 1x prefilled glass syringe of 1 or 1.5 mL volume, packed in steam sterilized blister.
- 1x 27G angled Cannula (S<sub>200</sub>-VISC 1%, 1.4%, 1.6%, 1.8%) sterilized by EO.
- 1x 25G angled Cannula (S<sub>200</sub>-VISC 2%, 3%) sterilized by EO.
- This Information for Use leaflet.
- Product identification labels for administrative purposes.

#### INDICATIONS

S<sub>200</sub>-VISC is indicated for use as a space occupying tissue protective fluid and lubricant during various intraocular surgical procedures including intra-/ extracapsular, cataract extraction, IOL implants and corneal transplantations. In surgical procedures in the anterior segment of the eye, instillation of the S<sub>200</sub>-VISC serves to maintain a deep anterior chamber during surgery, allowing for efficient manipulation with fewer traumas to the corneal endothelium and other surrounding tissues.

#### CONTRAINDICATIONS

S<sub>200</sub>-VISC must not be administered to patients who are known to be allergic to one of the constituents.

#### ADVERSE REACTIONS

The risks and side effects that may result from the use of S<sub>200</sub>-VISC are the same as those resulting from the use of any ophthalmic viscoelastic fluid. Sodium Hyaluronate is a physiological substance present in many tissues of the body. It is extremely well tolerated in human eyes. However, there have been reports of transient postoperative ocular inflammation (oral and/or topical steroid treatments were administered) and transient postoperative increase of intraocular pressure is likely to occur if the S<sub>200</sub>-VISC is not removed as completely as possible. Unexpected adverse reactions and complications related to S<sub>200</sub>-VISC must be reported to SIDAPHARM.

#### WARNINGS / PRECAUTIONS

- This product should only be used by qualified people who are familiar with the use of such a product.
- S<sub>200</sub>-VISC is intended only for intraocular surgery.
- It must not be re-sterilized.
- S<sub>200</sub>-VISC is intended for single use only.
- Do not administer after the expiry date.
- Opened and/or damaged sterile packs must not be used.
- To prevent the formation of air bubbles, remove the cap from the syringe by twisting it off carefully.
- Care should be taken to prevent the introduction of trapped air bubbles into the anterior chamber.
- Use only if solution is clear.

##### Interaction with other agents:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. S<sub>200</sub>-VISC should therefore not come into contact with operating instruments rinsed with these solutions or with ophthalmic preparations containing ammonium compounds as a preservative. This solution may appear cloudy or form precipitate when it comes into contact with quaternary ammonium compounds. If this is observed, the physician should remove the cloudy or precipitated material by irrigation and/or aspiration.

#### PREPARATION & ADMINISTRATION GUIDE

S<sub>200</sub>-VISC should be held at room temperature for approximately 30 minutes before use. The solution is administered through a narrow-gauge butt cannula. The amount required depends on the type of surgical procedure. There is a risk of the product drying out if it is left open for a long time during use.

- Remove the syringe from its packaging in a sterile environment
- Twist counterclockwise and remove the cap from the luer lock adapter of the syringe.
- Lock into position the enclosed cannula by twisting it clockwise
- Depress the plunger and discard the first 0.1 to 0.3 mL of fluid
- Do not overfill the eye chamber with S<sub>200</sub>-VISC.
- The product must be removed completely by irrigation and/or aspiration at the end of the operation. Otherwise, a mechanical blockage in the trabecular system could occasionally develop and cause a temporary postoperative rise in intraocular pressure. In patients diagnosed preoperatively with wide-angle glaucoma, pronounced myopia, diabetic retinopathy or uveitis, there is a high risk of an increase in intraocular pressure. If intraocular pressure rises to an excessive level, suitable treatment should be given to reduce the pressure.

#### STORAGE CONDITIONS

Store between 2-25 °C in a cool, dry place, and protect from light, heat and frost.

#### APPLICABLE REF.NUMBERS

Sodium Hyaluronate:	1.0%	1.4%	1.6%	1.8%	2.0%	3.0%
1.0 mL Ref.No.:	10000	10005	10006	10002	10009	10007
1.5 mL Ref.No.:	10010	10014	10016	10018	-	10030

#### DISCLAIMER

The manufacturer does not bear any responsibility for improper model selection by the surgeon, for improper handling, use, surgical technique applied or for any other iatrogenic error caused by the surgeon. This product is subject to change with or without prior notice. Improvement changes may be made in specification, shape and material. The manufacturer makes no expressed or implied warranty in connection with the sale of this product. The manufacturer is liable only for the design and production of this product and not for any incidental, consequential, indirect or exemplary damages of any kind, directly or indirectly arising from the purchase or use of it.

DE	DEUTSCH
----	---------

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