ENGLISH

PHACO-TIPS AND I/A-INSTRUMENTS (MULTIPLE USE)

INSTRUCTIONS FOR USE

These instructions for use are to be read and used only in conjunction with the manufacturer's instruction manual. Failure to consult the manufacturer's manual could result in improper use of this product.

Put pressure on the tip and screw the phaco tip onto the handpiece up to the stop position. Attention: Overtighten the tip can damage the multiple use / single use wrench as well as it can cause residue from the wrench on the surface of the tip.

TREATMENT ACCORDING TO EN ISO 17664

General principles: Before every use, all instruments must be cleaned, disinfected and sterilized; this applies particularly to the first use after delivery since all instruments are delivered unsterilized (cleaning and disinfection after removal of protective shipping packaging; sterilization). An effective cleaning and disinfection is an absolute necessity for effective sterilization. As part of your responsibility for the sterility of instruments during use, please ensure that only appropriate device- and product-specific procedures for cleaning/disinfection are used, that the devices used (disinfectant, sterilizer) are regularly maintained and inspected and that the va-

lided parameters are adhered to with every cycle. Only qualified personnel should perform the treatment. In addition, please observe the applicable legal regulations in your country as well as the hygiene regulations of the doctor's practice or hospital. This is especially true for the different guidelines regarding an effective inactivation of prions. **CLEANING AND DISINFECTION** **Principles according to EN ISO 15883:** For the cleaning and disinfection of instruments, a mechanical process (disinfectant) should be used, if possible. Due to its clearly reduced effectiveness, a manual procedure – even when using an ultrasound bath – should only be used if a mechanical process is not available. The use of a manual cleaning and disinfection procedure must be guaranteed by an additional product- and procedure-specific validation and is the user's responsibility. In both cases, pretreatment must be performed.

PRETREATMENT

Immediately after use (within no more than 1 hour), rough impurities must be removed from the products. To do so, use running water or a disinfectant solution; the disinfectant should be free of aldehyde (otherwise, possible fixing of contaminants), have a tested effectiveness (e.g., German Hygiene and Microbiology Society [DGHM], or FDA approval or CE marking), be suitable for instrument disinfection and be compatible with the instruments (see chapter "Material Durability"). For manual removal of impurities, only use a soft brush or a clean soft cloth that is only used for this purpose; never use metal brushes or steel wool. Please note that the disinfectant used in pretreatment is only for protecting persons and must not replace the disinfection step to be carried out later, after cleaning has been done. Check all instruments for visible contamination. Repeat the pretreatment if contamination is found.


MACHINE CLEANING/DESINFECTION

In selecting the disinfectant, make sure that • the disinfectant multiple use has a tested effectiveness (e.g. DGHM, or FDA approval or CE marking), • a tested program for thermal disinfection (at least 10 minutes at 93 °C) is used (with chemical disinfection there is the danger of disinfectant residue on the instruments), if possible, • the program used is suitable for the instruments and has sufficient rinsing cycles, • for washing down, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water (e.g. Aqua purificata) is used (EN 285, Appendix B), • the air used for drying is filtered and • the disinfectant is regularly maintained and inspected. **In selecting the cleaning agent used,** • it is basically suitable for cleaning the instruments (e.g. contains no abrasive materials), • if no thermal disinfection is used, a suitable disinfectant with a tested effectiveness (e.g. DGHM, or FDA approval

or CE marking) is also used and that it is compatible with the cleaning agent used and • the chemicals used are compatible with the instruments (see chapter "Material Durability"). The concentrations indicated by the manufacturer of the cleaning agent and, if necessary, of the disinfectant, must absolutely be complied with. **Sequence:** 1. Place the instruments in the disinfectant. While doing so, make sure that the instruments do not touch each other. 2. Start the program. 3. After the end of the program, remove the instruments from the disinfectant. 4. Inspect and pack the instruments as quickly as possible after removal (see "Inspection", "maintenance" and "packing" sections, if necessary after additional drying in a clean place).

Cleaning and disinfection device or other measures must ensure sufficient drying. **INSPECTION AND TESTING** Check all instruments after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping and contamination and sift out damaged instruments. Instruments that are still dirty must be re-cleaned and disinfected. **MAINTENANCE** Instrument oil or other oil may not be used.

PHACO-TIPS AND I/A-INSTRUMENTS (MULTIPLE USE)

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