





Tränenwegskanülen

Entwickelt zum Sondieren und Spülen der Tränenwege.

	<p>Tränenwegskanüle, gebogen Abgerundete stumpfe Spitze 1 VPE (10 Stück)</p> <table border="0"> <tr> <td>21 G x 28 mm (1¹/₈")</td> <td>Art.-Nr. 8000-21</td> <td rowspan="2">⊗</td> <td rowspan="2">STERILE</td> </tr> <tr> <td>21 G x 38 mm (1¹/₂")</td> <td>Art.-Nr. 8000-21 (38)</td> </tr> <tr> <td>23 G x 28 mm (1¹/₈")</td> <td>Art.-Nr. 8000-23</td> <td></td> <td></td> </tr> <tr> <td>25 G x 28 mm (1¹/₈")</td> <td>Art.-Nr. 8000</td> <td></td> <td></td> </tr> <tr> <td>25 G x 38 mm (1¹/₂")</td> <td>Art.-Nr. 8000 (38)</td> <td></td> <td></td> </tr> <tr> <td>26 G x 28 mm (1¹/₈")</td> <td>Art.-Nr. 8000-26</td> <td></td> <td></td> </tr> <tr> <td>26 G x 32 mm (1¹/₄")</td> <td>Art.-Nr. 8000-26 (32)</td> <td></td> <td></td> </tr> <tr> <td>26 G x 38 mm (1¹/₂")</td> <td>Art.-Nr. 8000-26 (38)</td> <td></td> <td></td> </tr> </table>	21 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-21	⊗	STERILE	21 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000-21 (38)	23 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-23			25 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000			25 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000 (38)			26 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-26			26 G x 32 mm (1 ¹ / ₄ ")	Art.-Nr. 8000-26 (32)			26 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000-26 (38)		
21 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-21	⊗	STERILE																												
21 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000-21 (38)																														
23 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-23																														
25 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000																														
25 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000 (38)																														
26 G x 28 mm (1 ¹ / ₈ ")	Art.-Nr. 8000-26																														
26 G x 32 mm (1 ¹ / ₄ ")	Art.-Nr. 8000-26 (32)																														
26 G x 38 mm (1 ¹ / ₂ ")	Art.-Nr. 8000-26 (38)																														
	<p>Tränenwegskanüle, gerade Abgerundete stumpfe Spitze 12 mm vom Ende von 19 G auf 23 G verjüngend 1 VPE (10 Stück)</p> <p>Art.-Nr. 8000S</p> <p>⊗ STERILE</p>																														

Hersteller für alle Produkte dieser Seite: Anodyne Surgical